

PriorRG: Prior-Guided Contrastive Pre-training and Coarse-to-Fine Decoding for Chest X-ray Report Generation


Kang Liu, Zhuoqi Ma, Zikang Fang, Yunan Li, Kun Xie, Qiguang Miao*

Xidian University, China

kangliu@stu.xidian.edu.cn, qgmiao@xidian.edu.cn



西安电子科技大学
XIDIAN UNIVERSITY




PriorRG
Code & Checkpoints



Awesome Radiology
Report Generation



Wechat




AAAI-26 / IAAI-26 / EAAI-26
JANUARY 20-27, 2026 | SINGAPORE

Introduction

- Background:** Chest X-ray report generation aims to reduce radiologists' workload by automatically producing high-quality preliminary reports.
- Contribution:** ① We propose **PriorRG**, which integrates **patient-specific prior knowledge** to generate **context-aware and disease progression-oriented reports**. ② We introduce a **prior-guided contrastive pre-training scheme** that mirrors diagnostic reasoning by using clinical context to guide spatiotemporal features extraction, thereby improving alignment with report semantics and boosting medical image-text retrieval performance. ③ We present a **prior-aware coarse-to-fine decoding** that incrementally integrates clinical context, disease progression patterns, and hierarchical visual cues, enhancing the clinical accuracy and fluency of generated reports.
- Results:** Extensive experiments on MIMIC-CXR and MIMIC-ABN datasets demonstrate that PriorRG outperforms state-of-the-art methods, achieving a **3.6% BLEU-4** and **3.8% F1 score** improvement on MIMIC-CXR, and a **5.9% BLEU-1** gain on MIMIC-ABN. Furthermore, PriorRG is computationally efficient, with an average inference time of **0.1245 seconds per sample** and a GPU memory footprint of **7.68 GB**.

Motivation

PriorRG jointly models **disease progression (DP)** and **clinical context (CC)** to enhance **cross-modal alignment** and improve the quality of generated reports.




Report Generation Performance

Dataset	Method	Venue	NLG Metrics ↑						CE Metrics ↑		
			B-1	B-2	B-3	B-4	MTR	R-L	P	R	F1
M-CXR	KiUT	CVPR'23	0.393	0.243	0.159	0.113	0.160	0.285	0.371	0.318	0.321
	METransformer	CVPR'23	0.386	0.250	0.169	0.124	0.152	0.291	0.364	0.309	0.311
	CoFE	ECCV'24	-	-	-	0.125	0.176	0.304	0.489	0.370	0.405
	DCG	MM'24	0.397	0.258	0.166	0.126	0.162	0.295	0.441	0.414	0.404
	MAN	AAAI'24	0.396	0.244	0.162	0.115	0.151	0.274	0.411	0.398	0.389
	R2GenGPT	Meta-Radio'23	0.411	0.267	0.186	0.134	0.160	0.297	0.392	0.387	0.389
	Med-LLM	MM'24	-	-	-	0.128	0.161	0.289	0.412	0.373	0.395
	R2-LLM	AAAI'24	0.402	0.262	0.180	0.128	0.175	0.291	0.465	<u>0.482</u>	<u>0.473</u>
	SEI	MICCAI'24	0.382	0.247	0.177	0.135	0.158	0.299	<u>0.523</u>	0.410	0.460
	HERGen	ECCV'24	0.395	0.248	0.169	0.122	0.156	0.285	-	-	-
M-ABN	MPO	AAAI'25	0.416	<u>0.269</u>	<u>0.191</u>	<u>0.139</u>	0.162	<u>0.309</u>	0.436	0.376	0.353
	PriorRG (Ours)	-	<u>0.412</u>	0.290	0.220	0.175	0.189	0.324	0.541	0.485	0.511
	△(%) ↑	-	-0.4	+2.1	+2.9	+3.6	+1.3	+1.5	+1.8	+0.3	+3.8
	CMN [◇]	ACL'21	0.256	0.147	0.095	0.066	0.110	0.230	0.466	0.454	0.460
	SEI [◇]	MICCAI'24	<u>0.267</u>	<u>0.157</u>	<u>0.104</u>	<u>0.073</u>	<u>0.114</u>	<u>0.231</u>	0.466	0.408	0.435
	PriorRG (Ours)	-	0.326	0.201	0.139	0.102	0.140	0.242	0.467	0.476	0.471
	△(%) ↑	-	+5.9	+4.4	+3.5	+2.9	+2.6	+1.1	+0.1	+2.2	+1.1

Table 1: Comparison of report generation performance with SOTA methods on MIMIC-CXR (M-CXR) and MIMIC-ABN (M-ABN) datasets. Δ indicates the performance difference between PriorRG and the best baseline. $\diamond$ denotes reproduced results; others are cited from the original papers. The best and second-best values are in bold and underlined, respectively.

Medical Image-text Retrieval Results




Ablation Study of Report Generation on MIMIC-CXR


Model	Stage 1	Stage 2			NLG Metrics ↑						CE Metrics ↑		
		CC	PI	Hidden States	B-1	B-2	B-3	B-4	MTR	R-L	P	R	F1
(a)	✓	×	×	×	0.355	0.219	0.149	0.108	0.154	0.262	0.521	0.431	0.472
(b)	✓	×	×	✓	0.357	0.216	0.143	0.102	0.154	0.255	0.513	0.480	0.496
(c)	✓	✓	×	×	0.405	0.283	0.214	0.170	0.186	0.320	0.517	0.460	0.487
(d)	✓	✓	✓	×	0.404	0.285	0.217	0.173	0.187	0.328	0.540	0.477	0.507
(e)	✓	✓	✓	✓	0.400	0.282	0.214	0.171	0.186	0.323	0.539	0.465	0.499
(f)	×	✓	✓	✓	0.387	0.271	0.206	0.165	0.180	0.320	0.526	0.408	0.459
PriorRG	✓	✓	✓	✓	0.412	0.290	0.220	0.175	0.189	0.324	0.541	0.485	0.511

Table 2: Ablation study on MIMIC-CXR dataset. “CC” and “PI” denote clinical context and prior image, respectively.

Overview



Expert Alignment & Generalization Capability



Effect of Prior-Knowledge Ratio and Progressive Fusion

Prior Type	Ratio	NLG Metrics ↑		CE Metrics ↑		
		B-2↑	F1↑	P	R	F1
Prior Image (PI)	0%	0.287	0.513			
	10%	0.287	0.513			
	30%	0.288	0.513			
	50%	0.289	0.515			
	100%	0.289	0.515			
Clinical Context (CC)	0%	0.139	0.492			
	10%	0.157	0.495			
	30%	0.187	0.498			
	50%	0.220	0.501			
	100%	0.294	0.507			

Table 9: Report generation performance on the MIMIC-CXR dataset under varying proportions of Prior Image (PI) and Clinical Context (CC).

Qualitative Analysis

Input Data	Radiologist-generated Report	AI-generated Report by SEI	AI-generated Report by PriorRG (Ours)
Current Image Prior Image	Comparison is made to the prior study from [1]. The feeding tube, left IJ catheter and endotracheal tube are unchanged in position. There is persistent cardiomegaly. There is unchanged left retrocardiac opacity. There are no signs for overt pulmonary edema. There is a small right-sided pleural effusion as well. Overall, these findings are stable.	in comparison with the earlier study of this date, there has been placement of an endotracheal tube with its tip approximately 5 cm above the carina. otherwise, little change.	Comparison is made to prior study from [1]. The endotracheal tube, feeding tube, and left IJ central line are unchanged in position. There is unchanged cardiomegaly. There is a persistent left retrocardiac opacity and left-sided pleural effusion which is stable. There are no signs for overt pulmonary edema. No pneumothoraces are identified.
Current Image Prior Image	As compared to the previous radiograph, the Dobhoff tube shows now normal course. The tip projects over the middle parts of the stomach. No complications, notably no pneumothorax. Otherwise, the image is unchanged.	as compared to the previous radiograph, the dobhoff tube has been repositioned, the course of the tube is unremarkable, the tip of the tube is not included on the image. otherwise, there is no relevant change.	As compared to the previous radiograph, the patient has received a Dobhoff catheter. The course of the catheter is unremarkable, the tip of the catheter projects over the middle parts of the stomach. There is no evidence of complications, notably no pneumothorax. Otherwise, the radiograph is unchanged.